

# Nitrocellulose Analysis by Low Angle Light Scattering

Paul Clarke, Viscotek Europe

NC Meeting November 2001

# Topics

- Background
- Brief review of Light Scattering theory
- LALS detector
- NC data by LALS compared to RALS method
- Discussion/Conclusions
- Further work

## Principles of Light Scattering Measurements

Light Scattering from Solutions is governed by the Rayleigh Eqn.

$$[1] \quad \frac{KC}{R_0} = \frac{1}{Mw} + 2A_2C$$

Where  $R_0$  is the Rayleigh Ratio at zero scattering angle,  
 $M_w$  is the weight-average Molecular Weight,  
 $C$  is the concentration of the solution,  
 $K$  is an optical constant,  
 $A_2$  is the 2nd virial coefficient.


Note that at low concentration(GPC conditions?) the 2nd term on the right becomes negligible.

$$[1a] \quad R_0 = KCMw$$


## Principles of LS: Rayleigh Ratio

Scattered light cannot be measured at zero angle, of course, so it must be obtained from estimations based on measurement(s) at some higher angle  $\theta$ .  $R_\theta$  is defined as the ratio of scattered light intensity per unit volume of the solution to the incident light intensity, normalized.

$$[2 - 1] \quad R_\theta = \frac{I_\theta r^2}{I_0 V}$$




# Principles of Light Scattering - Angular Dependence



Small Molecule (radius <15nm): Light is scattered equally at all angles.

Zero angle only true result



Large Molecule (radius >15nm): Light is scattered less at higher angles.

For larger macromolecules the phenomenon of *intramolecular interference* produces ~~asymmetry~~ **anisotropy** in the scattered light. Destructive interference occurs at backward scattering angles.

© Viscotek 2001

## Principles of LS: Angular Dependence

If the dimensions of the scattering molecules are much smaller than the wavelength of the incident light, the molecule behaves virtually as a point source of scattered light. The result is a symmetrical scattering pattern and  $R_\theta$  may be set equal to  $R_0$  without significant error.

$$R_\theta \approx R_0$$

The accuracy of this equation is dependent on  $\theta$  and on the size of the macromolecule. For  $\theta = 90^\circ$ , for example, this is a good approximation for  $R_g < 15$  nm, which corresponds roughly to Mw of 150,000 for chain polymers.

Note that for spherical molecules such as proteins,  $R_g$  of 15 nm would correspond to Mw of several million!

## Principles of LS: Angular Correction

For large macromolecules exhibiting LS disymmetry, one must estimate  $R_0$  from measurement(s) of  $R_\theta$ . We can define a *particle scattering function*  $P_\theta$  to be

$$[8 - 1] \quad P_\theta = \frac{R_\theta}{R_0}$$

There are three ways to obtain  $R_0$  or, equivalently,  $P_\theta$ . These ways correlate with the commercial light scattering instruments.


1. Measure  $R_\theta$  at very low angles. - LALS
2. Measure  $R_\theta$  at several angles and extrapolate. - MALS
3. Measure  $R_\theta$  at  $90^\circ$  and calculate  $P_\theta$  using viscosity data.

- RALS/Viscometer  
 VISCOTEK  
New Technology in GPC and Viscosities


© Viscotek 2001

# Angular Dependence

## *Stiffer Chain Molecules ( $a = 0.8$ )*



# LALLS DETECTOR - CUTAWAY





 New Technology in GPC and Viscosities

 **Viscotek**  
New Technology in Viscosities

 VISCOTEK 2001

# LALS-RALS


## *Response for Small and Large Molecules*




# Experimental

- Calibrate system with Narrow PS
- Check calibration with broad PS standard
- Run 4 nitrocellulose solutions supplied by Alan Macdonald from AWE
- Calculate data by both LALS and RALS/Visc


# PS Standard by RALS & LALS




# AWE Sample 04




# AWE Sample 08



# AWE Sample 09



# AWE Sample 07



## Analysis Summary (Triple Detector)

File	Sample	Mw	Mn	PD	IVw	Rg (nm)	alpha	log K	dn/dc	Conc	% Rec.	Cal
4-2	PS 250K Broad - RALS	254,500	96,720	2.63	0.8138	18.17	0.717	-3.933	0.185	1.510	0.000	nc_rals
4-2L	PS 250K Broad - LALS	250,300	95,030	2.63	0.8127	18.07	0.716	-3.923	0.185	1.510	0.000	nc_lals
5-1	AWE 04	200,600	56,170	3.57	3.542	26.49	0.914	-4.282	0.085	2.028	0.000	nc_rals
5-1L	AWE 04	208,800	67,130	3.11	3.555	27.01	0.959	-4.545	0.085	2.028	0.000	nc_lals
6-1	AWE 08	238,500	112,000	2.13	8.235	38.12	0.839	-3.570	0.121	1.764	0.000	nc_rals
6-1L	AWE 08	246,200	137,000	1.80	8.255	38.98	0.945	-4.172	0.121	1.764	0.000	nc_lals
7-1	AWE 09	443,800	102,200	4.34	5.881	39.01	0.809	-3.748	0.088	0.760	0.000	nc_rals
7-1L	AWE 09	417,600	141,400	2.95	6.136	39.57	0.952	-4.551	0.088	0.760	0.000	nc_lals
8-1	AWE 07	334,300	95,910	3.49	3.793	32.46	0.864	-4.169	0.073	1.576	0.000	nc_rals
8-1L	AWE 07	328,500	118,200	2.78	3.825	32.59	0.954	-4.674	0.073	1.576	0.000	nc_lals

GPC TDA300 - Triple mode

column GMHx1 2x


solvent THF

**Viscotek Europe Ltd**


*Triple Detector GPC analysis*

*Jul 12, 2001*


# AWE08, LALS data




# Mark-Houwink Plot



## Mark-Houwink Plot



# Light Scattering Chromatograms



# Conclusions

- LALS data broadly confirms RALS data.
- Some small differences in molecular weight but mainly <5%
- LALS raises further questions about ‘pre-peak’

# Further work

- More structured study:
  - Analysis of ‘fresh’ solutions
  - Analysis with time
- Focus on ‘pre-peak’
  - Measure MW separately by LALS
  - $d\eta/dc$ ?

# Thanks!

- Nicola Marheineke
- Alan Macdonald
- Max Haney